

KMTY

Gourmet Tepan Yaki

for a unique cooking experience on Miele induction cooktops.

EAN: 4002516349846 / Material number: 11555900 / Old Material Number: 26996100EU1



- For preparing light and healthy dishes
- With FixPosition – four silicone feet ensure a secure grip
- Even frying thanks to premium multi-layer stainless steel plate.
- Ideal for cooktops with PowerFlex and full-surface induction
- Not suitable for gas and electric cooktops
- External dimensions: WxDxH 9 1/2"x16"x7/8" /frying surface: WxD 8 1/4"x13 3/8"