

TXR860WP Eco & Steam



EAN: 4002516379836 / Material number: 11620030 /
Old Material Number: 12XR8602USA

Design	
Stackable	•
Side-by-Side	•
Slot-in	•
Door hinging	right, Reversible
Convertible door hinge	•
Type	
Product brand	T1 White Edition
Heat-pump dryer	•
Design	
Appliance color	Lotus white
Control panel color	Lotus white
Door design	Lotus white with chrome ring
Fascia panel version	5° inclined
Control type	Touch control
Display	MTouch
Drying results	
PerfectDry precision drying	•
FragranceDos ² fragrance system	•
Intelligent drum reversal	•
Gentle laundry care	
Honeycomb drum	•
Load in lbs. (kg)	8.0
User convenience	
Networking with Miele@home	•
Delay start up to 24 hours	•
Countdown indicator	•
Time display	•
Programme sequence indicator	Yes
Drying assistant	•
Information text	•
Favorite programs	•
Drum lighting	LED
MultiLingua	•
Condensed water container	•
Integrated condensed water drainage	•
AddLoad	•
Efficiency and sustainability	
EcoFeedback	•
EcoDry technology	•
Maintenance-free heat exchanger	•
ProfiEco motor	•
Drying programs	
Delicates	•
Express	•
Outerwear	•
Proofing	•

TXR860WP Eco & Steam



EAN: 4002516379836 / Material number: 11620030 /
Old Material Number: 12XR8602USA

Normal	•
Wrinkle-free	•
Smoothing	•
Woolens/Silks	•
Dress Shirts	•
Time Dry	•
Denims	•
Table Linens/Drapes	•
Bed Linens	•
Cold Air	•
Basket	•
Sportswear	•
Soft Steam	•
Pillows	•
Drying extras	
Anti-crease action	•
Gentle tumble	•
Freshen up	•
Acoustic signal	•
Quality	
Soft ribs	•
Safety	
Lock code	•
Empty out container" indicator	•
Clean filter" indication	•
Optical interface	•
Technical data	
Dimensions (width) in in. (mm)	23 1/2 (596)
Dimensions (height) in in. (mm)	33 1/2 (850)
Dimensions (depth) in in. (mm)	25 3/8 (643)
Depth with door open in in. (mm)	42 1/2 (1077)
Weight in lbs. (kg)	135 (61)
Total connected load in kW	1.00
Voltage in V	120
Fuse rating in A	15
Frequency in Hz	60
Length of supply lead in ft. (m)	6 (2)
Length of supply lead in ft. (m)	6 (2.0)
Length of supply lead in ft. (m)	6 (2)
Type of coolant	R134a
Refrigerant quantity in lbs (kg)	1.03 (0.47)
Standard accessories	
Fragrance flacon	•
Voucher for a second fragrance flacon	•

TXR860WP Eco & Steam



EAN: 4002516379836 / Material number: 11620030 /
Old Material Number: 12XR8602USA

Available display languages

Available display languages via MultiLingua

bahasa malaysia, dansk,
deutsch, english (AU),
english (CA), english (GB),
english (US), español,
français (CA), français
(FR), hrvatski, italiano,
magyar, nederlands, norsk,
polski, português, русский,
română, slovenčina,
slovenščina, srpski, suomi,
svenska, türkçe, українська,
čeština, ελληνικά,
[?][?][?][?][?][?], [?][?], [?]
[?][?], [?][?][?][?], [?][?][?]