


APTD 500

Basket

for optimum loading of up to 8 bottles and other hollow containers.



- Easy to repurpose using different baskets – simple conversion
- Hygienic reprocessing thanks to individual injector nozzles
- Supplied with adapter
- Max. load clearance 14 3/4" (375 mm), diameter 4" (100 mm)

EAN: 4002516381709 / Material number: 11622900 / Old Material Number: 67650200D

Product category	
Tank dishwashers	•
Machine affiliation	
PTD 702	•
PTD 703	•
PTD 704	•
Product category	
Basket for flasks and other hollow containers	•
Product characteristics	
Material	Wire, rilsan-coated
Color	White
Capacity	
Flasks up to 14 3/4" (375 mm) high [number]	8
Dimensions and weight	
External dimensions, net height in inches	8 15/16 (227)
External dimensions, net width in inches	19 11/16 (500)
External dimensions, net depth in inches	19 11/16 (500)
External dimensions, gross height in inches	18 3/4 (475)
External dimensions, gross width in inches	21 5/16 (540)
External dimensions, gross depth in inches	21 5/16 (540)
Net weight in lbs	7 (2.93)
Gross weight in lbs	10 (4.40)