

A 814/1

Holder

for optimum reprocessing of tips for air scaler and piezo scaler.



- For use in A 315/1 module, as well as A 105/1, O 177/1, and APWD 061 upper
- Independent of manufacturer, for a number of different tips
- Use of several holders possible in combination with A 813/1
- Use in thermal disinfectors PWD 8531 and PWD 8532 possible

EAN: 4002516428558 / Material number: 11691870 / Old Material Number: 68681401D

Product category	
Thermal disinfectors	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8528	•
G 7831	•
PG 8536	•
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
PWA 8672	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Holder	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of dental instruments	•
Capacity	
Ultrasonic tips [number]	1
Dimensions and weight	
External dimensions, net height in inches	1 15/16 (48)
External dimensions, net width in inches	0 7/8 (22)
External dimensions, gross height in inches	4 5/16 (108)
External dimensions, gross width in inches	4 5/8 (117)
External dimensions, gross depth in inches	12 1/16 (305)
Net weight in lbs	2 (0.76)
Gross weight in lbs	3 (0.99)