

A 865
Holder
for perfectly holding multi-functional and 3F syringes.



- For use in A 105/1 upper basket and A 315/1 module with A 800 tubular filter
- Compatible with syringes from several well-known manufacturers

EAN: 4002516413752 / Material number: 11709860 / Old Material Number: 69686500D

Product category	
Thermal disinfectors	•
Cleaning machine, small, dental	•
Large-capacity washer-disinfectors, dental	•
Machine affiliation	
PG 8581	•
PG 8591	•
PWA 8672	•
PWD 7121	•
PWD 7122	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Holder	•
Application area	
Reprocessing of dental instruments	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Dimensions and weight	
External dimensions, net height in inches	4 13/16 (121)
External dimensions, net width in inches	1 1/8 (28)
External dimensions, gross height in inches	4 5/16 (108)
External dimensions, gross width in inches	4 5/8 (117)
External dimensions, gross depth in inches	12 1/16 (305)
Net weight in lbs	1 (0.30)
Gross weight in lbs	2 (0.94)