

A 865  
Holder  
for perfectly holding multi-functional and 3F syringes.



- For use in A 105/1 upper basket and A 315 module
- Compatible with syringes from several well-known manufacturers

EAN: 4002516413752 / Material number: 11709860 / Old Material Number: 69686500D

<b>Product category</b>	
Thermal disinfectors	•
<b>Machine affiliation</b>	
PG 8581	•
PG 8591	•
<b>Product category</b>	
Holder	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
<b>Application area</b>	
Reprocessing of dental instruments	•
<b>Dimensions and weight</b>	
External dimensions, net height in inches	4 13/16 (121) in (mm)
External dimensions, net width in inches	1 1/8 (28) in (mm)
External dimensions, net depth in inches	0 (0) in (mm)
External dimensions, gross height in inches	4 5/16 (108)
External dimensions, gross width in inches	4 5/8 (117)
External dimensions, gross depth in inches	12 1/16 (305)
Net weight in lbs	1 (0.30)
Gross weight in lbs	3 (0.93)