

ProCare Shine 12 GC - 80 Tabs
80 detergent tabs
for optimum cleaning of dishes, cutlery and glasses.

EAN: 4002516420590 / Material number: 11712580 / Old Material Number: 61914201USA



- Particularly simple handling thanks to water-soluble wrapping
- Free from phosphates and chlorine
- For light to medium soiling
- For use in commercial dishwashers

Product category	
Fresh-water dishwashers	•
ProfiLine dishwashers	•
Applications	
Removal of coffee and tea residues	Good
Removal of starch/protein	Good
Removal of dried food residues	Good
Removal of burnt-on grease	Good
Removal of grease	Very good
Dispensing recommendation	
Dispensing recommendation per rinse	1 Tab
Intensity of soiling	Light to moderate
pH value of solution	10.2
Recommended water range in mmol/l	0-0.53
Product characteristics	
Medium type	Mildly alkaline Cleaning agent
Aggregate form	Tab
Water-soluble film	•
Suitable for glasses	•
Suitable for cutlery	•
Suitable for plastic	•
Suitable for stainless steel	•
Suitable for porcelain/ceramics	•
Suitable for aluminum	•
Suitable for silver	•
Gentle on decorative material	•
Protects against corrosion	•
Prevents limescale and protects wash items	•
Does not contain phosphates	•
Dimensions and weight	
External dimensions, net height in inches	9 1/16 (230)
External dimensions, net width in inches	7 1/2 (190)
External dimensions, net depth in inches	3 11/16 (93)
External dimensions, gross height in inches	9 1/16 (230)
External dimensions, gross width in inches	3 9/16 (90)
External dimensions, gross depth in inches	7 1/2 (190)
Number of detergent tabs	80
Load size in g	45.15 (1280)
Gross weight in lbs	0.04 (1)
Compliance, test certificates, and symbols	
GHS/CLP compliant	•