

A 830/1

Mesh Basket

ideal holder for separable mouth specula.



- For use in insert E 337/1 and E 802/1
- 3 compart. Ø 8 mm, 14 compart. 10 x 19 mm, 2 compart. 8 x 19 mm
- Mouth specula handles can be placed in compartments of inserts
- Suitable for mouth specula ≥ 3

EAN: 4002516429203 / Material number: 11729480 / Old Material Number: 68683001D

<b>Product category</b>	
Thermal disinfectors	•
<b>Machine affiliation</b>	
PG 8528	•
G 7831	•
PG 8536	•
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
<b>Product category</b>	
Insert/module for dental instruments	•
<b>Product characteristics</b>	
Material	Plastic
Color	Black
Number of compartments	16
<b>Application area</b>	
Reprocessing of dental instruments	•
<b>Dimensions and weight</b>	
External dimensions, net height in inches	0 9/16 (14) in (mm)
External dimensions, net width in inches	1 15/16 (49) in (mm)
External dimensions, net depth in inches	6 5/16 (159) in (mm)
External dimensions, gross height in inches	2 11/16 (67)
External dimensions, gross width in inches	6 3/4 (170)
External dimensions, gross depth in inches	9 1/16 (230)
Net weight in lbs	1 (0.02)
Gross weight in lbs	1 (0.02)