

Boost CX1 Cat & Dog
 Bagless canister vacuum cleaners
 for the highest hygiene demands, in a compact design.



- Powerful cleaning performance thanks to Vortex technology
- Always stays on track thanks to TrackDrive
- Maximum air hygiene with HEPA filter: *Ideal for allergy sufferers
- with additional TurboBrush – ideal for pet lovers
- Convenient height adjustment thanks to Comfort Telescopic Tube

* According to EN 60312-1

EAN: 4002516420019 / Material number: 11735550 /
 Old Material Number: 41NCE031USA

General information	
Model identifier	Boost CX1 Cat & Dog
Manufacturer material number	11735550
EAN number	4002516420019
Design	
Bagless vacuum cleaner	•
Product line	Boost CX1
Motor type	PowerLine
Product type	SNCEO
Appliance type and design	
Appliance colour	Obsidian black
Cleaning	
Vortex Technology	•
User convenience	
Electronic suction power setting	4-stage rotary selector
Ergonomic carrying handle	•
Ergonomic handle	•
Telescopic suction wand	Comfort telescopic wand
Park system for storing	on one side
Mobility	
TrackDrive	•
Operating radius in ft (m)	32 (10)
Cord length in ft (m)	21 (6.5)
Filter system	
Exhaust filter	HEPA AirClean filter
Dustbag volume in cu. ft. (l)	0.04 (1.0)
Technical data	
Power	1200
Volume in dB	78.0
Width in inches (mm)	15 3/4 (400)
Height in inch (mm)	11 (280)
Depth in inches (mm)	11 (280)
Weight in lbs (kg)	18 (7.8)
Standard accessories	
Upholstery nozzle	•
Crevice nozzle	•
Dust brush	•
Floorhead SBD 365-3	•
Turbobrush STB 305-3 (DOUBLE; LOCKED)	•
Cleaning brush	No