

A 813/1

Adapter

optimum reprocessing of transmission instruments with an external spray channel.



- For use in A 105/1 upper baskets and in APWD 061 baskets
- Two-piece, comprising a top and bottom part incl. accessories

EAN: 4002516432555 / Material number: 11763880 / Old Material Number: 68681301D

Product category	
Thermal disinfectors	•
Machine affiliation	
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
Product category	
Holder and/or adapter for transmission instruments	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of dental instruments	•
Capacity	
Transmission instruments [number]	4
Dimensions and weight	
External dimensions, net height in inches	120
External dimensions, net width in inches	120
External dimensions, net depth in inches	300
External dimensions, gross height in inches	4 5/16 (108)
External dimensions, gross width in inches	4 5/8 (117)
External dimensions, gross depth in inches	12 1/16 (305)
Net weight in lbs	1 (0.19)
Gross weight in lbs	3 (1.00)
Companion pack	
E 362 blind stopper [number]	3
Nozzle with silicone hose [number]	4