


APLW 097

Lower basket

For loading various inserts for wide-neck laboratory glassware.

EAN: 4002516429807 / Material number: 11764830 / Old Material Number: 69309700D



- Useable surface: 490 x 470 mm (19 5/16 X 18 1/2 in)
- Vertical clearance with upper basket: 180 mm (7 3/32 in)
- Vertical clearance without upper basket: 460 mm (18 1/8 in)

<b>Machine/accessory categorization</b>	
PLW 8505	•
<b>Product category</b>	
Lower basket	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, gross height in inches	7 7/8 (200)
External dimensions, gross width in inches	21 3/4 (550)
External dimensions, gross depth in inches	23 5/8 (600)
Net weight in lbs	9 (3.71)
Gross weight in lbs	13 (5.52)