


APLW 098

Upper basket

With 36 connections for injector nozzles.



- Injector nozzles are not included with the machine
- Vertical clearance: 160 mm (6 5/16 In)
- Suitable for use with lower basket APLW 095 and APLW 097

EAN: 4002516430216 / Material number: 11764840 / Old Material Number: 69309800D

| | |
|---|--------------------------|
| Machine/accessory categorization | |
| PLW 8505 | • |
| Product category | |
| Upper basket | • |
| Product characteristics | |
| Material | Stainless steel, Plastic |
| Color | Blue, Stainless Steel |
| Application area | |
| Reprocessing of laboratory glassware | • |
| Capacity | |
| Narrow-neck glassware [number] | 36 |
| Dimensions and weight | |
| External dimensions, gross height in inches | 7 7/8 (200) |
| External dimensions, gross width in inches | 21 3/4 (550) |
| External dimensions, gross depth in inches | 23 5/8 (600) |
| Net weight in lbs | 10 (4.32) |
| Gross weight in lbs | 14 (6.12) |