

DKFS 31-R

Reactivatable Longlife AirClean charcoal filter
for Miele DAS 2x20, DAS 4x20, and DAS 4x30 Recirculation Hoods.



EAN: 4002516441892 / Material number: 11793540 /
Old Material Number: 28996436EU1

- Set for initial commissioning
- Air so fresh and odor-free that you can feel it in your kitchen
- Can easily be regenerated in the oven every six months
- Service life 3 years
- Special construction type for low flow resistance
- Contains: 2 DKF 31-R filters with frames