

CK/1 DIN 51

Canister key

for easy opening of ProCare 5 l, 10 l, and 20 l canisters.

EAN: 4002516480839 / Material number: 11851600 / Old Material Number: 66990002EU1



Product category	
ProCare	•
Dimensions and weight	
External dimensions, net height in inches	1 1/8 (27)
External dimensions, net width in inches	8 13/16 (223)
External dimensions, net depth in inches	3 3/8 (85)
External dimensions, gross height in inches	2 3/16 (55)
External dimensions, gross width in inches	7 1/8 (180)
External dimensions, gross depth in inches	9 7/8 (250)
Gross weight in lbs	0 (0)

- Sturdy and durable manual key
- Easy to open and close media canisters without force
- Can be used on both sides