

APLW 215

Injector module

for the optimum loading of large-volume plastic bottles up to 50 L.



- For use on level 1
- Equipped with a special injector nozzle
- Equipped with a holding frame
- Max. bottle diameter: 1 3/16" (380 mm)
- Bottle height: 22 1/16" to 24 7/16" (560 to 620 mm)
- Bottle neck diameter: min. 1 1/4" (30 mm) internal, max. 3 1/8" (80 mm) external

EAN: 4002516525981 / Material number: 11853940 / Old Material Number: 69621500D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 7111	•
Product category	
Other inserts/other modules	•
Module for large-volume laboratory glassware	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Large-volume laboratory glassware [number]	1
Dimensions and weight	
External dimensions, net height in inches	22 1/4 (565)
External dimensions, net width in inches	18 1/2 (469)
External dimensions, net depth in inches	20 3/4 (526)
External dimensions, gross height in inches	23 1/4 (590)
External dimensions, gross width in inches	20 1/2 (520)
External dimensions, gross depth in inches	20 7/8 (530)
Net weight in lbs	9 (3.84)
Gross weight in lbs	12 (5.14)