


APLW 106

Basket

for the optimum loading of injector modules or inserts.



- For use on levels 1 to 4
- Basket with two docking connections
- For loading up to 2 injector modules or inserts
- Automatic closing of the connections when not in use
- Clearances dependent on which baskets are used and their position

EAN: 4002516524717 / Material number: 11854610 / Old Material Number: 69610600D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 7111	•
Product category	
Basket with module connection	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Part of the SimpleLoad system	•
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	6 (151)
External dimensions, net width in inches	22 1/16 (559)
External dimensions, net depth in inches	24 1/2 (621)
External dimensions, gross height in inches	6 15/16 (175)
External dimensions, gross width in inches	22 1/2 (570)
External dimensions, gross depth in inches	24 13/16 (630)
Net weight in lbs	7 (3.05)
Gross weight in lbs	9 (3.75)