

WA UC 1501 L

UltraColor liquid detergent 1.5 l

highly effective liquid detergent for colors and dark/black items.

EAN: 4002516531593 / Material number: 11979450 / Old Material Number: 11997231USA



- Excellent spot removal even at low temperatures
- Contains 7 enzymes for ultimate cleanliness
- Bottles made from 100% recycled materials – eco-friendly
- ColorProtect for long-lasting color intensity
- Free from micro-plastics and biodegradable
- Miele Aqua: a refreshingly natural fragrance