

WA SO 1503 L
UltraSoft fabric softener 1.5 l
for soft, bouncy laundry and a fresh Aqua fragrance.

EAN: 4002516531470 / Material number: 11987230 / Old Material Number: 11997233USA



- Natural fragrance of pure cleanliness and freshness
- Minimizes creasing and makes ironing laundry easier
- Bottles made from 100% recycled materials – eco-friendly
- Free from micro-plastics and biodegradable
- Prevents electrostatic charges in washer dryers
- Very economical – for 50 wash loads