

APLW DS7111

Retrofitting kit

liquid dispensing for retrofitting of the PLW 7111 with a third dispenser pump.

EAN: 4002516648697 / Material number: 11993930 / Old Material Number: 69242001D



- Can be retrofitted by Customer Service
- For the use of additional process chemicals
- A maximum of two retrofitting kits can be installed

| Product category | |
|--|-------------------------------|
| SlimLine laboratory glasswashers | • |
| Machine affiliation | |
| PLW 7111 | • |
| PLW 7111 [Safety] | • |
| PWD 7121 | • |
| PWD 7122 | • |
| Product category | |
| Conversion kit | • |
| Dispenser system | • |
| Product characteristics | |
| Material | Stainless steel, Plastic |
| Color | Stainless Steel, Green, Black |
| Application area | |
| Reprocessing of laboratory glassware | • |
| Reprocessing of dental instruments | • |
| Reprocessing of lumened instruments | • |
| Reprocessing of trays | • |
| Reprocessing of air regulators | • |
| Reprocessing of breathing masks | • |
| Reprocessing of compressed-air unit carrier frames | • |
| Reprocessing of compressed air bottles | • |
| Dimensions and weight | |
| External dimensions, gross height in inches | 3 3/4 (95) |
| External dimensions, gross width in inches | 6 3/4 (170) |
| External dimensions, gross depth in inches | 21 11/16 (550) |
| Net weight in lbs | 1 (0.42) |
| Gross weight in lbs | 2 (0.64) |