

FA C 452 L

Set of 3 x Cocoon fragrance flacon
for a pleasant fragrance – soothing and restful.

EAN: 4002516561446 / Material number: 12023770 / Old Material Number: 11997229USA



- Individually adjustable fragrance intensity at any time
- For a fresh fragrance that lasts up to four weeks
- Beautifully fragrant laundry – without fabric softener
- Specially developed for Miele by perfumers from in the Provence
- For Miele tumble dryers and the laundry cabinet with FragranceDos
- For 150 drying cycles or 210 applications in the laundry cabinet