

FA DF 452 L
Set of 3x Miele DryFresh™
for 150 drying cycles – Freshen-up and fragrance.

EAN: 4002516561422 / Material number: 12026080 / Old Material Number: 11997236USA



- With Freshplex™ to neutralize unpleasant odors
- Freshens up dry laundry in the “Timed Drying” program
- Beautifully fragrant laundry – without fabric softener
- Can be combined with other Miele Fragrance Flacons
- Specially developed for Miele by perfumers from in the Provence
- For all Miele T1 Tumble Dryers with FragranceDos