

APWD 325

Mesh Basket

for reprocessing hoses used in dental nitrous oxide sedation systems.



- Simultaneous reprocessing of up to 5 double hose sets possible
- With adapter set for secure connection to injector system
- Reliable internal/external cleaning and drying of hoses
- Option of parallel patient mask reprocessing in mesh tray E 146
- Large capacity maintained for other instruments in each cycle

EAN: 4002516610908 / Material number: 12036420 /
Old Material Number: 68632500D

Product category

Thermal disinfectors

•

Machine affiliation

PG 8581

•

PG 8591

•

PWA 8672

•

PWD 8682

•

PWD 8682 CD

•

PWD 8692

•

Product category

Other inserts

•

Application area

Reprocessing of dental instruments

•

Reprocessing of different utensils

•

Product characteristics

Material

Stainless steel, Silicone

Color

Stainless Steel, Grey

Dimensions and weight

External dimensions, net height in inches

5 1/4 (132)

External dimensions, net width in inches

12 13/16 (324)

External dimensions, net depth in inches

14 15/16 (379)

External dimensions, gross height in inches

6 5/16 (160)

External dimensions, gross width in inches

13 1/4 (335)

External dimensions, gross depth in inches

15 (380)

Net weight in lbs

5 (2.23)

Gross weight in lbs

6 (2.62)