


APWD 329

Mesh Basket

for reprocessing dental saliva extractors.



- Suitable for up to 80 dental saliva extractors
- 80 compartments, 3/4" x 3/4" (18 x 18 mm) each
- Mesh width on base 5/16" x 5/16" (8 x 8 mm)

EAN: 4002516596080 / Material number: 12036470 /
Old Material Number: 69632900D

Product category

Thermal disinfectors •
Cleaning machine, small, dental •

Machine affiliation

PG 8581 •
PG 8591 •
PWA 8672 •
PWD 8682 •
PWD 8682 CD •
PWD 8692 •

Product category

Insert/module for dental instruments •

Application area

Reprocessing of dental instruments •

Product characteristics

Material Stainless steel
Color Stainless steel

Capacity

Siphons [number] 80

Dimensions and weight

External dimensions, net height in inches 8 13/16 (223)
External dimensions, net width in inches 8 13/16 (223)
External dimensions, net depth in inches 3 3/16 (80)
External dimensions, gross height in inches 3 15/16 (100)
External dimensions, gross width in inches 9 1/2 (240)
External dimensions, gross depth in inches 9 1/2 (240)
Net weight in lbs 1 (0.58)
Gross weight in lbs 2 (0.71)