


DKF 31-900 P

Active AirClean charcoal filter

DAS 4x40, DAS 8x30, and DKFS 31-P for DAS 2x20, DAS 4x20, and DAS 4x30.

EAN: 4002516606277 / Material number: 12123400 / Old Material Number: 28997257USA



- Air so fresh and odor-free that you can feel it in your kitchen
- Effective neutralization of unpleasant cooking odors
- Pleated material reduces noise levels
- Quick and easy filter replacement
- Replace filter every 6 months to maintain performance
- Contains: 2 filters without frames