


DKFS 31-900 P

Active AirClean charcoal filter

for Miele DAS 2x20, DAS 4x20, and DAS 4x30 Recirculation Hoods.



EAN: 4002516606468 / Material number: 12123420 /
Old Material Number: 28997259USA

- Set for initial commissioning
- Air so fresh and odor-free that you can feel it in your kitchen
- Pleated material reduces noise levels
- Quick and easy filter replacement
- Replace filter every 6 months to maintain performance
- Contains: 2 DKF 31-P filters with frames