

APLW 870

Set of injector nozzles

with support, for cleaning sedimentation vials (Imhoff cones).



- Set consists of 4 special nozzles with holders
- Max. diameter 4 3/4" max (120 mm, max.) height 18 1/2" (470 mm)
- For use with separately available injector module A 300/2
- For use with separately available injector module A 300/3
- For use with separately available injector module A 620
- Part of the SimpleLoad system

EAN: 4002516638124 / Material number: 12151110 / Old Material Number: 69687000D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Part of the SimpleLoad system	•
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Measuring cylinders [number]	4
Dimensions and weight	
External dimensions, net height in inches	3 3/8 (85)
External dimensions, net width in inches	6 15/16 (175)
External dimensions, net depth in inches	11 7/16 (290)
External dimensions, gross height in inches	3 15/16 (100)
External dimensions, gross width in inches	7 1/8 (180)
External dimensions, gross depth in inches	11 13/16 (300)
Net weight in lbs	1 (0.39)
Gross weight in lbs	2 (0.56)