

APLW 870

Set of injector nozzles
with support, for cleaning sedimentation vials (Imhoff cones).



- Set consists of 4 special nozzles with holders
- Max. diameter 4 3/4" max (120 mm, max.) height 18 1/2" (470 mm)
- For use with separately available injector module A 300/2
- For use with separately available injector module A 300/3
- For use with separately available injector module A 280
- Part of the SimpleLoad system

EAN: 4002516638124 / Material number: 12151110 / Old Material Number: 69687000D

| | |
|---|--------------------------|
| Product category | |
| Washer-disinfectors, laboratory | • |
| Large-capacity washer-disinfectors, laboratory | • |
| SlimLine laboratory glasswashers | • |
| Machines for cleaning industrial parts | • |
| Large-chamber laboratory glasswashers for cleaning industrial parts | • |
| SlimLine laboratory glasswashers for cleaning industrial parts | • |
| Machine affiliation | |
| PLW 8636 [IPC] | • |
| PLW 8636 [LAB / LAB MON] | • |
| PG 8504 | • |
| PG 8583 | • |
| PG 8583 CD | • |
| PG 8593 | • |
| PLW 8604 | • |
| PLW 8615 | • |
| PLW 8616 | • |
| PLW 8617 | • |
| PLW 8683 | • |
| PLW 8683 CD | • |
| PLW 8693 | • |
| PLW 7111 | • |
| PLW 7111 [IPC] | • |
| Product category | |
| Injector nozzle | • |
| Application area | |
| Reprocessing of laboratory glassware | • |
| Product characteristics | |
| Material | Stainless steel, Plastic |
| Color | Stainless Steel, Black |
| Part of the SimpleLoad system | • |
| Capacity | |
| Measuring cylinders [number] | 4 |
| Dimensions and weight | |
| External dimensions, net height in inches | 3 3/8 (85) |
| External dimensions, net width in inches | 6 15/16 (175) |
| External dimensions, net depth in inches | 11 7/16 (290) |
| External dimensions, gross height in inches | 3 15/16 (100) |
| External dimensions, gross width in inches | 7 1/8 (180) |
| External dimensions, gross depth in inches | 11 13/16 (300) |
| Net weight in lbs | 1 (0.40) |
| Gross weight in lbs | 1 (0.56) |