

ProCare Dent 10 MA - 5 l

Liquid detergent, mildly alkaline, 5 l
for optimum reprocessing of dental instruments.



- Specially designed for cleaning dental instruments
- For use with Dental Washer-Disinfectors
- Very gentle on material
- Liquid concentrate

EAN: 4029577143428 / Material number: 12186860 / Old Material Number: 66913203EU1

Product category	
Thermal disinfectors	•
Applications	
Removal of blood	Good
Removal of organic residues	Good
Dispensing recommendation	
Dispensing recommendation per rinse	3-5 ml/l (0.3-0.5%)
pH value of solution	10.6-10.7
Main applications	
Dental practices	•
Product characteristics	
Medium type	Mildly alkaline Cleaning agent
Aggregate form	Liquid concentrate
Ingredients	Enzymes, Anionic tensides, Non-ionic tensides
Suitable for dental handpieces	•
Very gentle on materials	•
Suitable for standard instruments and transmission instruments	•
Dimensions and weight	
External dimensions, net height in inches	8 9/16 (217)
External dimensions, net width in inches	5 5/8 (142)
External dimensions, net depth in inches	9 13/16 (248)
External dimensions, gross height in inches	9 7/8 (250)
External dimensions, gross width in inches	7 1/2 (190)
External dimensions, gross depth in inches	5 15/16 (150)
Gross weight in lbs	0.18 (5)
Fill volume in ml	169.07 (5000) fl oz (US) (Milliliter)
Compliance, certificates/symb.	
CE	•
GHS/CLP compliant	•