

ProCare Dent 11 A - 10 kg
 Powder cleaning agent, alkaline, 10 kg
 for optimum reprocessing of standard instruments.



- Suitable for dental instruments
- For use with G 7881 Dental Washer-Disinfector
- High cleaning performance
- Powder

EAN: 4029577143442 / Material number: 12186990 / Old Material Number: 66911202EU1

| | |
|--|-------------------------------|
| Product category | |
| Thermal disinfectors | • |
| Applications | |
| Removal of blood | Good |
| Removal of organic residues | Good |
| Dispensing recommendation | |
| Dispensing recommendation per rinse | 3 g/l (30 g/cycle) |
| pH value of solution | 11.1 |
| Main applications | |
| Dental practices | • |
| Product characteristics | |
| Medium type | Alkaline cleaning agent |
| Aggregate form | Powder |
| Ingredients | Phosphate, Non-ionic tensides |
| Good cleaning capability | • |
| High degree of material protection against corrosion | • |
| Suitable for standard instruments | • |
| Dimensions and weight | |
| External dimensions, net height in inches | 10 9/16 (268) |
| External dimensions, net width in inches | 10 7/16 (265) |
| External dimensions, net depth in inches | 10 7/16 (265) |
| External dimensions, gross height in inches | 11 1/16 (280) |
| External dimensions, gross width in inches | 10 11/16 (270) |
| External dimensions, gross depth in inches | 10 11/16 (270) |
| Load size in g | 352.73 (10000) |
| Gross weight in lbs | 0.35 (10) |
| Compliance, certificates/symb. | |
| CE | • |
| GHS/CLP compliant | • |