

ProCare Dent 30 C - 5 l

Neutralizing agent, acidic, 5 l

for optimum neutralization of alkaline residues, based on organic acids.



- Suitable for dental instruments
- For use with Dental Washer-Disinfectors
- Very gentle on material
- Liquid concentrate

EAN: 4029577143466 / Material number: 12187250 / Old Material Number: 66915302EU1

<b>Product category</b>	
Thermal disinfectors	•
<b>Applications</b>	
Removal of limescale and other acid-soluble impurities	Good
<b>Dispensing recommendation</b>	
Dispensing recommendation per rinse	1 ml/l (0.1%)
pH value of solution	3.0
<b>Main applications</b>	
Dental practices	•
<b>Product characteristics</b>	
Medium type	Acidic neutralising agent
Aggregate form	Liquid concentrate
Ingredients	Citric acid
Does not contain surfactants	•
Very good material compatibility	•
Suitable for standard instruments and transmission instruments	•
<b>Dimensions and weight</b>	
External dimensions, net height in inches	8 9/16 (217)
External dimensions, net width in inches	5 5/8 (142)
External dimensions, net depth in inches	9 13/16 (248)
External dimensions, gross height in inches	9 7/8 (250)
External dimensions, gross width in inches	7 1/2 (190)
External dimensions, gross depth in inches	5 15/16 (150)
Gross weight in lbs	0.21 (6)
Fill volume in ml	169.07 (5000) fl oz (US) (Milliliter)
<b>Compliance, certificates/symb.</b>	
CE	•
GHS/CLP compliant	•