

ProCare Dent 30 C - 5 l
Neutralizing agent, acidic, 5 l
for optimum neutralization of alkaline residues, based on organic acids.



- Suitable for dental instruments
- For use with Dental Washer-Disinfectors
- Very gentle on material
- Liquid concentrate

EAN: 4029577143466 / Material number: 12187250 / Old Material Number: 66915302EU1

| | |
|--|--|
| Product category | |
| Thermal disinfectors | • |
| Applications | |
| Removal of limescale and other acid-soluble impurities | Very good |
| Dispensing recommendation | |
| Dispensing recommendation per rinse | 1 ml/l (0.1%) |
| Dispensing recommendation (per rinse) | 1 ml/l (0.1%) |
| Large opening for siphon | • |
| pH value at 68°F (20°C) | 1,0 |
| Suitable for all levels of water hardness | • |
| Main applications | |
| Dental practices | • |
| Product characteristics | |
| Medium type | Acidic neutralising agent |
| Aggregate form | Liquid concentrate |
| Ingredients | Citric acid, Preservatives |
| Minimum storage temperature in °F/°C | -3 |
| Maximum storage temperature in °F/°C | 30 |
| Dimensions and weight | |
| External dimensions, gross depth in inches (mm) | 150.00 |
| Product characteristics | |
| Medium type | Acidic neutralising agent |
| Aggregate form | Liquid concentrate |
| Ingredients | Citric acid, Preservatives |
| Minimum storage temperature in °F/°C | -3 |
| Maximum storage temperature in °F/°C | 30 |
| Does not contain surfactants | • |
| Suitable for sensitive materials | • |
| Neutralization of alkaline carryover | • |
| Excellent rinsing properties | • |
| Prevention of limescale formation in wash cabinet | • |
| Does not contain phosphorus compounds | • |
| Suitable for standard instruments and transmission instruments | • |
| Suitable for anodized aluminum | • |
| Very good material compatibility | • |
| Dimensions and weight | |
| External dimensions, net width in inches | 5 5/8 (142) |
| External dimensions, net depth in inches | 9 13/16 (248) |
| External dimensions, gross height in inches | 9 7/8 (250) |
| External dimensions, gross width in inches | 7 1/2 (190) |
| External dimensions, gross depth in inches | 5 15/16 (150) |
| Gross weight in lbs | 0.21 (6) |
| Fill volume in ml | 169.07 (5000) fl oz (US) (Milliliter) |
| Compliance, test certificates, and symbols | |
| CE | • |
| GHS/CLP compliant | • |