

ProCare Dent 40 - 1 l [Typ 1]

Rinse aid, 1 l

with pouring aid, for final rinsing of the wash items.



- Improves and speeds up drying of dental instruments
- For use with Dental Washer-Disinfectors
- Easier dispensing thanks to precise pouring aid
- Liquid concentrate

EAN: 4029577155117 / Material number: 12187410 / Old Material Number: 66917102EU1

Product category	
Thermal disinfectors	•
Dispensing recommendation	
Dispensing recommendation per rinse	0.3 ml/l (0.03%)
Dispensing recommendation (per rinse)	0.3 ml/l (0.03%)
pH value at 68°F (20°C)	2,0
Suitable for all levels of water hardness	•
Small opening as pouring aid	•
Main applications	
Dental practices	•
Product characteristics	
Medium type	Rinsing agent
Aggregate form	Liquid concentrate
Ingredients	Preservatives, Surfactants, Citric acid
Minimum storage temperature in °F/°C	0
Maximum storage temperature in °F/°C	30
Dimensions and weight	
External dimensions, gross depth in inches (mm)	105.00
Product characteristics	
Medium type	Rinsing agent
Aggregate form	Liquid concentrate
Ingredients	Preservatives, Surfactants, Citric acid
Minimum storage temperature in °F/°C	0
Maximum storage temperature in °F/°C	30
Suitable for porcelain/ceramics	•
Improves and shortens the drying process	•
Neutralization of alkaline carryover	•
Suitable for standard instruments and transmission instruments	•
Suitable for all metals used in medical devices	•
Dimensions and weight	
External dimensions, net height in inches	9 1/2 (241)
External dimensions, net width in inches	4 3/16 (106)
External dimensions, net depth in inches	2 5/16 (58)
External dimensions, gross height in inches	9 11/16 (245)
External dimensions, gross width in inches	2 3/8 (60)
External dimensions, gross depth in inches	4 3/16 (105)
Gross weight in lbs	0.04 (1)
Fill volume in ml	33.81 (1000) fl oz (US) (Milliliter)
Compliance, test certificates, and symbols	
CE	•
GHS/CLP compliant	•