


APTD 221
Basket set
for mixed wash items in bistros and cafés.



- Base basket incl. 3 inserts for plates, glasses, and cutlery
- Highly versatile
- Base basket dimensions: 19 11/16" x 19 11/16" (500 x 500 mm)

EAN: 4002516717812 / Material number: 12305290 / Old Material Number: 67022101D

Product category	
Tank dishwashers	•
Throughfeed tank dishwashers	•
Machine affiliation	
PTD 901	•
PFD 400	•
PFD 400 U	•
PFD 401	•
PFD 401 U	•
PFD 402	•
PFD 402 U	•
PFD 404	•
PFD 404 U	•
PFD 405	•
PFD 405 U	•
PFD 407	•
PFD 407 U	•
PFD 408 U	•
PTD 702	•
PTD 703	•
PTD 704	•
Product category	
Basket for mixed tableware	•
Product characteristics	
Material	Plastic
Color	Grey
Number of rows	3
Dimensions and weight	
External dimensions, net height in inches	0 (0)
External dimensions, net width in inches	0 (0)
External dimensions, net depth in inches	0 (0)
External dimensions, gross height in inches	8 7/8 (225)
External dimensions, gross width in inches	21 11/16 (550)
External dimensions, gross depth in inches	21 11/16 (550)
Net weight in lbs	5 (1.92)
Gross weight in lbs	7 (2.96)