

APLW 875
Retrofitting kit
for retrofitting a DI water booster pump to a PLW 7111.



- Installed by a service technician

EAN: 4002516792444 / Material number: 12385630 /
Old Material Number: 69242005D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 7111	•
PLW 7111 [Safety]	•
Product category	
Other accessories	•
Application area	
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Grey
Dimensions and weight	
External dimensions, gross height in inches	4 7/8 (123)
External dimensions, gross width in inches	12 7/16 (315)
External dimensions, gross depth in inches	19 15/16 (505)
Net weight in lbs	5 (2.13)
Gross weight in lbs	5 (2.46)