

ProCare Universal 81 C - 250 g
Machine descaler, acidic, 250 g
for removing persistent limescale deposits.



- Removal of limescale and contamination trapped within
- Prolongs service life of heater elements
- Specially developed for Miele Professional machines

EAN: 4002516770862 / Material number: 12456130 / Old Material Number: 66920601V1

Product category	
Fresh water dishwashers	•
Fresh-water dishwashers	•
ProfiLine dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Washing machines	•
Laboratory glasswashers	•
Large-chamber laboratory glasswashers	•
Thermal disinfectors	•
Large-chamber washer-disinfectors	•
Washer-disinfectors, dental	•
Laboratory glasswasher, medical	•
Large-capacity washer-disinfectors, medical	•
Applications	
Removal of limescale and other acid-soluble impurities	•
Dispensing	
Recommended dosage per wash cycle	1–2 doses
Suitable for all levels of water hardness	•
Manual dispensing	•
Product characteristics	
Medium type	Machine descaler
Aggregate form	Powder
Ingredients	Citric acid
pH value at 68°F (20°C)	4,5
Does not contain phosphorus compounds	•
Very gentle on materials	•
Very good material compatibility	•
Free from microplastics	•
Storage information	
Minimum storage temperature in °F/°C	5
Maximum storage temperature in °F/°C	40
Dimensions and weight	
External dimensions, net height in inches	4 15/16 (124)
External dimensions, net width in inches	2 5/8 (66)
External dimensions, net depth in inches	2 5/8 (66)
External dimensions, gross height in inches (mm)	4 3/4 (120)
External dimensions, gross width in inches (mm)	2 13/16 (70)
External dimensions, gross depth in inches (mm)	2 13/16 (70)
Load size in g	8.82 (250)
Gross weight in lbs (kg)	1 (0.33)
Compliance, test certificates, and symbols	
GHS/CLP compliant	•