

APLW 895
Retrofitting kit
drain water cooling for retrofitting in the PLW 7111.



- For cooling drain water in temperature-sensitive drain
- Installed by Miele Service

EAN: 4002516888062 / Material number: 12616020 / Old Material Number: 69242016D

Product category	
SlimLine laboratory glasswashers	•
Machine affiliation	
PLW 7111	•
Product category	
Conversion kit	•
Application area	
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel, Plastic
Dimensions and weight	
External dimensions, gross height in inches	11 7/16 (290)
External dimensions, gross width in inches	4 1/16 (102)
External dimensions, gross depth in inches	12 13/16 (324)
Net weight in lbs	2 (0.92)
Gross weight in lbs	2 (1.09)