

A 100

Upper basket

For the optimum loading of injector modules or inserts.



- Upper basket with two module connections
- for holding up to 2 injector modules or inserts
- Automatic closing of the connections when not in use
- Loading height with module with plastic supports 7.5 inches
- Clearances: width 19 1/4" (489 mm), length 17 5/16" (440 mm)

EAN: 4002516841296 / Material number: 12642640 / Old Material Number: 69610002D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Machines for cleaning industrial parts	•
<b>Machine affiliation</b>	
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
<b>Product category</b>	
Basket with module connection	•
Upper basket	•
<b>Application area</b>	
Reprocessing of pipettes/vials	•
Reprocessing of laboratory glassware	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
<b>Dimensions and weight</b>	
External dimensions, net height in inches	5 9/16 (141)
External dimensions, net width in inches	20 11/16 (525)
External dimensions, net depth in inches	20 13/16 (528)
External dimensions, gross height in inches	6 3/4 (170)
External dimensions, gross width in inches	21 7/8 (555)
External dimensions, gross depth in inches	22 1/4 (565)
Net weight in lbs	6 (2.92)
Gross weight in lbs	8 (3.56)