


A 105/1

Upper basket

Height adjustable, without injector nozzles, for various instruments.



- Right-hand side: space for 22 injectors
- Left-hand side free for inserts
- #Active internal drying of handpieces and contra-angle handpieces
- Clearances H 225 +/- 30, W 360, D 445 mm
- A 800 tubular filter can be used in both injector manifolds

EAN: 4002516854814 / Material number: 12664820 / Old Material Number: 69610502D

| Product category                            |                          |
|---------------------------------------------|--------------------------|
| Thermal disinfectors                        | •                        |
| Cleaning machine, small, dental             | •                        |
| Washer-disinfectors, medical                | •                        |
| Machine affiliation                         |                          |
| PG 8581                                     | •                        |
| PG 8562                                     | •                        |
| PG 8582                                     | •                        |
| PG 8582 CD                                  | •                        |
| PG 8591                                     | •                        |
| PG 8592                                     | •                        |
| PWA 8672                                    | •                        |
| PWD 8682                                    | •                        |
| PWD 8682 CD                                 | •                        |
| PWD 8692                                    | •                        |
| Product category                            |                          |
| Injector basket                             | •                        |
| Upper basket                                | •                        |
| Application area                            |                          |
| Reprocessing of ENT instruments             | •                        |
| Reprocessing of dental instruments          | •                        |
| Product characteristics                     |                          |
| Material                                    | Stainless steel, Plastic |
| Color                                       | Stainless Steel, Black   |
| Option to use a central filter              | •                        |
| Capacity                                    |                          |
| Transmission instruments [number]           | 22                       |
| Lumened instruments [number]                | 22                       |
| Dimensions and weight                       |                          |
| External dimensions, net height in inches   | 8 (202)                  |
| External dimensions, net width in inches    | 20 13/16 (528)           |
| External dimensions, net depth in inches    | 20 13/16 (528)           |
| External dimensions, gross height in inches | 8 11/16 (220)            |
| External dimensions, gross width in inches  | 21 5/16 (540)            |
| External dimensions, gross depth in inches  | 21 5/16 (540)            |
| Net weight in lbs                           | 11 (4.88)                |
| Gross weight in lbs                         | 12 (5.57)                |
| Accessory pack                              |                          |
| A 838 key set [number]                      | 1                        |