

A 861 Holding frames for A 301/4 and A 301/5 injector modules.



- Arrangement: 3 x 6 fields
- For use on A 301/4 and A 301/5 injector modules
- For the optimum positioning of laboratory glassware
- Part of the SimpleLoad system

EAN: 4002516854876 / Material number: 12665440 / Old Material Number: 69686101D

Product category

Washer-disinfectors, laboratory	•
SlimLine laboratory glasswashers	•
Machines for cleaning industrial parts	•
SlimLine laboratory glasswashers for cleaning industrial parts	•

Machine affiliation

PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
PLW 7111 [IPC]	•

Product category

Other inserts/other modules	•
-----------------------------	---

Application area

Reprocessing of laboratory glassware	•
--------------------------------------	---

Product characteristics

Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Number of compartments	18
Compartment width in inches	2 15/16 (74)
Compartment length in inches	2 15/16 (74)
Part of the SimpleLoad system	•

Dimensions and weight

External dimensions, net height in inches	3 3/16 (80)
External dimensions, net width in inches	7 7/8 (200)
External dimensions, net depth in inches	13 13/16 (350)
External dimensions, gross height in inches	5 1/8 (130)
External dimensions, gross width in inches	10 1/16 (255)
External dimensions, gross depth in inches	19 1/8 (485)
Net weight in lbs	2 (0.59)
Gross weight in lbs	3 (0.94)