

A 202

Mobile unit

for the optimum loading of 4 DIN mesh trays.



- Suitable for 4 DIN mesh trays on 2 levels
- Insert frames: Removable for the reprocessing of tall wash items
- 6 additional connection options via Luer lock connection
- With integrated spray arm
- Loading volumes for the different levels differ from each other

EAN: 4002516854043 / Material number: 12667130 / Old Material Number: 69620201D

| Product category                            |                          |
|---------------------------------------------|--------------------------|
| Laboratory glasswashers                     | •                        |
| Washer-disinfectors, medical                | •                        |
| Machines for cleaning industrial parts      | •                        |
| Machine affiliation                         |                          |
| PLW 8636 [IPC]                              | •                        |
| PLW 8636 [LAB / LAB MON]                    | •                        |
| PG 8562                                     | •                        |
| PG 8582                                     | •                        |
| PG 8582 CD                                  | •                        |
| PG 8583                                     | •                        |
| PG 8583 CD                                  | •                        |
| PG 8592                                     | •                        |
| PG 8593                                     | •                        |
| PLW 8683                                    | •                        |
| PLW 8683 CD                                 | •                        |
| PLW 8693                                    | •                        |
| PWD 8682                                    | •                        |
| PWD 8682 CD                                 | •                        |
| PWD 8692                                    | •                        |
| Product category                            |                          |
| Mobile unit for inserts                     | •                        |
| Mobile unit for mesh trays                  | •                        |
| Application area                            |                          |
| Reprocessing of laboratory glassware        | •                        |
| Reprocessing of surgical instruments        | •                        |
| Product characteristics                     |                          |
| Material                                    | Stainless steel, Plastic |
| Color                                       | Stainless Steel, Black   |
| Capacity                                    |                          |
| DIN mesh trays [number]                     | 4                        |
| Dimensions and weight                       |                          |
| External dimensions, net height in inches   | 8 13/16 (223)            |
| External dimensions, net width in inches    | 20 7/8 (529)             |
| External dimensions, net depth in inches    | 21 3/8 (542)             |
| External dimensions, gross height in inches | 10 11/16 (270)           |
| External dimensions, gross width in inches  | 22 1/2 (570)             |
| External dimensions, gross depth in inches  | 21 11/16 (550)           |
| Net weight in lbs                           | 13 (6.05)                |
| Gross weight in lbs                         | 15 (6.76)                |