


E 103/1

Mesh Basket

For the optimum loading of test tubes up to 12 x 75 mm.



- Suitable for about 200 Test tubes (up to 12 x 75 mm)
- Divided into 6 sections
- The mesh width is 8 x 8 mm
- Including lid A 13

EAN: 4002516854104 / Material number: 12667340 / Old Material Number: 69510303D


| Product category                                                    |                 |
|---------------------------------------------------------------------|-----------------|
| Laboratory glasswashers                                             | •               |
| Large-chamber laboratory glasswashers                               | •               |
| SlimLine laboratory glasswashers                                    | •               |
| Machines for cleaning industrial parts                              | •               |
| Large-chamber laboratory glasswashers for cleaning industrial parts | •               |
| SlimLine laboratory glasswashers for cleaning industrial parts      | •               |
| Machine affiliation                                                 |                 |
| G 7825                                                              | •               |
| PG 8527                                                             | •               |
| PG 8535                                                             | •               |
| PG 8536                                                             | •               |
| PLW 8636 [IPC]                                                      | •               |
| PLW 8636 [LAB / LAB MON]                                            | •               |
| PG 8504                                                             | •               |
| PG 8583                                                             | •               |
| PG 8583 CD                                                          | •               |
| PG 8593                                                             | •               |
| PLW 8615                                                            | •               |
| PLW 8616                                                            | •               |
| PLW 8617                                                            | •               |
| PLW 8683                                                            | •               |
| PLW 8683 CD                                                         | •               |
| PLW 8693                                                            | •               |
| PLW 6011                                                            | •               |
| PLW 6111                                                            | •               |
| PLW 7111                                                            | •               |
| PLW 7111 [IPC]                                                      | •               |
| PLW 8505                                                            | •               |
| Product category                                                    |                 |
| Insert/module for test tubes                                        | •               |
| Application area                                                    |                 |
| Reprocessing of test tubes                                          | •               |
| Reprocessing of laboratory glassware                                | •               |
| Product characteristics                                             |                 |
| Material                                                            | Stainless steel |
| Color                                                               | Stainless steel |
| Capacity                                                            |                 |
| Test tubes [number]                                                 | 200             |
| Dimensions and weight                                               |                 |
| External dimensions, net height in inches                           | 4 13/16 (122)   |
| External dimensions, net width in inches                            | 7 7/8 (200)     |
| External dimensions, net depth in inches                            | 8 11/16 (220)   |
| External dimensions, gross height in inches                         | 5 9/16 (140)    |
| External dimensions, gross width in inches                          | 9 11/16 (245)   |
| External dimensions, gross depth in inches                          | 9 7/8 (250)     |

E 103/1

Mesh Basket

For the optimum loading of test tubes up to 12 x 75 mm.

EAN: 4002516854104 / Material number: 12667340 / Old Material Number: 69510303D



|                   |          |
|-------------------|----------|
| Net weight in lbs | 1 (0.56) |
|-------------------|----------|

|                     |          |
|---------------------|----------|
| Gross weight in lbs | 2 (0.73) |
|---------------------|----------|

|                       |  |
|-----------------------|--|
| <b>Accessory pack</b> |  |
|-----------------------|--|

|                   |  |
|-------------------|--|
| A 13 lid [number] |  |
|-------------------|--|

|  |   |
|--|---|
|  | 1 |
|--|---|