


E 104/1

Mesh Basket

for the optimum loading of test tubes up to 12 x 105 mm.



- Suitable for about 200 Test tubes (up to 12 x 105 mm)
- Capacity including lid A 13
- The mesh width is 8 x 8 mm

EAN: 4002516854111 / Material number: 12667370 / Old Material Number: 69510403D


Product category	
Laboratory glasswashers	•
Large-chamber laboratory glasswashers	•
SlimLine laboratory glasswashers	•
Machines for cleaning industrial parts	•
Large-chamber laboratory glasswashers for cleaning industrial parts	•
SlimLine laboratory glasswashers for cleaning industrial parts	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 7111 [IPC]	•
PLW 8505	•
Product category	
Insert/module for test tubes	•
Application area	
Reprocessing of test tubes	•
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel
Color	Stainless steel
Capacity	
Test tubes [number]	200
Dimensions and weight	
External dimensions, net height in inches	6 (152)
External dimensions, net width in inches	7 7/8 (200)
External dimensions, net depth in inches	9 1/16 (230)
External dimensions, gross height in inches	6 1/2 (165)
External dimensions, gross width in inches	9 1/2 (240)
External dimensions, gross depth in inches	9 7/8 (250)

E 104/1

Mesh Basket

for the optimum loading of test tubes up to 12 x 105 mm.

EAN: 4002516854111 / Material number: 12667370 / Old Material Number: 69510403D



Net weight in lbs	1 (0.62)
-------------------	----------

Gross weight in lbs	2 (0.84)
---------------------	----------

Accessory pack	
-----------------------	--

A 13 lid [number]	1
-------------------	---