

A 300/2

Injector module

For optimum loading of laboratory glassware – can be equipped with 8 nozzles.



- Suitable for lab. glassware such as Erlenmeyer and round flasks.
- Module w/o injector nozzles – equipped based on individual needs
- For up to 8 injector nozzles
- Part of the SimpleLoad system
- Nozzle spacing 4 5/8" across the width, 4 5/8" across the depth
- Particularly suitable for A 840 and A 841 (nozzle diameter 1/4" (6 mm))

EAN: 4002516853954 / Material number: 12668780 /
Old Material Number: 69630004D

Product category	
Laboratory glasswashers	•
SlimLine laboratory glasswashers	•
Machines for cleaning industrial parts	•
SlimLine laboratory glasswashers for cleaning industrial parts	•
Machine affiliation	
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
PLW 7111 [IPC]	•
Product category	
Insert/module for narrow-neck glassware	•
Insert for various bottles	•
Application area	
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel
Color	Stainless steel, Black
Part of the SimpleLoad system	•
Capacity	
Narrow-neck glassware [number]	8
Dimensions and weight	
External dimensions, net height in inches	2 7/8 (72)
External dimensions, net width in inches	5 1/4 (133)
External dimensions, net depth in inches	18 7/8 (479)
External dimensions, gross height in inches	7 7/8 (200)
External dimensions, gross width in inches	20 11/16 (525)
External dimensions, gross depth in inches	21 5/16 (540)
Net weight in lbs	1 (0.45)
Gross weight in lbs	2 (0.84)