

A 301/4

Injector module

For optimum loading of laboratory glassware – can be equipped with 18 nozzles.

EAN: 4002516854739 / Material number: 12668790 / Old Material Number: 69630106D



- Suitable for lab. glassware such as Erlenmeyer and round flasks.
- Module w/o injector nozzles – equipped based on individual needs
- For up to 18 injector nozzles
- Part of the SimpleLoad system
- Nozzle spacing 3" (78 mm) across the width, 3" (78 mm) across the depth
- Particularly suitable for A 842 and A 843 (nozzle diameter 4 mm)

Product category	
Laboratory glasswashers	•
SlimLine laboratory glasswashers	•
Machines for cleaning industrial parts	•
SlimLine laboratory glasswashers for cleaning industrial parts	•
Machine affiliation	
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
PLW 7111 [IPC]	•
Product category	
Insert/module for narrow-neck glassware	•
Insert for various bottles	•
Application area	
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless steel, Black
Part of the SimpleLoad system	•
Capacity	
Narrow-neck glassware [number]	18
Dimensions and weight	
External dimensions, net height in inches	2 7/8 (72)
External dimensions, net width in inches	6 13/16 (173)
External dimensions, net depth in inches	18 7/8 (479)
External dimensions, gross height in inches	3 15/16 (100)
External dimensions, gross width in inches	8 1/8 (205)
External dimensions, gross depth in inches	19 3/16 (487)
Net weight in lbs	2 (0.90)
Gross weight in lbs	4 (1.89)