


A 302/2

Injector module

For optimum loading of laboratory glassware — can be equipped with 32 nozzles.



- Suitable for lab. glassware such as Erlenmeyer and round flasks.
- Module w/o injector nozzles — equipped based on individual needs
- For up to 32 injector nozzles
- Part of the SimpleLoad system
- Nozzle spacing 2 3/8" across the width, 2 1/4" across the depth
- Particularly suitable for A 844 and A 845 (nozzle diameter 2.5 mm)

EAN: 4002516854098 / Material number: 12669210 /
Old Material Number: 69630204D

Product category	
Laboratory glasswashers	•
SlimLine laboratory glasswashers	•
Machines for cleaning industrial parts	•
SlimLine laboratory glasswashers for cleaning industrial parts	•
Machine affiliation	
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
PLW 7111 [IPC]	•
Product category	
Insert/module for narrow-neck glassware	•
Insert for various bottles	•
Application area	
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless steel, Black
Part of the SimpleLoad system	•
Capacity	
Narrow-neck glassware [number]	32
Dimensions and weight	
External dimensions, net height in inches	2 7/8 (72)
External dimensions, net width in inches	7 11/16 (195)
External dimensions, net depth in inches	18 7/8 (479)
External dimensions, gross height in inches	3 9/16 (90)
External dimensions, gross width in inches	8 1/8 (205)
External dimensions, gross depth in inches	19 11/16 (500)
Net weight in lbs	2 (1.12)
Gross weight in lbs	4 (1.96)