

A 304
Module
for the optimum loading of up to 98 vials.



- Suitable for 98 centrifuge tubes, phials, test tubes, for example
- The loading heights differ depending on the positioning
- 1 unit can be inserted into each of the baskets A 100, A 150, and APLW 106

EAN: 4002516854982 / Material number: 12669230 /
Old Material Number: 69630401D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machines for cleaning industrial parts	•
SlimLine laboratory glasswashers for cleaning industrial parts	•
Machine affiliation	
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
PLW 7111 [IPC]	•
Product category	
Insert/module for test tubes	•
Application area	
Reprocessing of pipettes/vials	•
Reprocessing of test tubes	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Capacity	
Phials [number]	98
Dimensions and weight	
External dimensions, net height in inches	5 1/8 (130)
External dimensions, net width in inches	8 3/4 (222)
External dimensions, net depth in inches	18 9/16 (471)
External dimensions, gross height in inches	5 15/16 (150)
External dimensions, gross width in inches	9 1/16 (230)
External dimensions, gross depth in inches	19 5/16 (490)
Net weight in lbs	5 (2.30)
Gross weight in lbs	7 (3.15)