


A 601

Mesh Basket

For the addition of another level in mobile unit A 501.



- Without spray arm

EAN: 4002516854081 / Material number: 12669810 / Old Material Number: 69660102D

Product category

Large-capacity washer-disinfectors, laboratory

•

Machine affiliation

PLW 8615

•

PLW 8616

•

PLW 8617

•

Product category

Other inserts/other modules

•

Application area

Reprocessing of laboratory glassware

•

Product characteristics

Material

Stainless steel

Color

Stainless Steel

Capacity

DIN mesh trays [number]

3

SPRI mesh trays [number]

2

Dimensions and weight

External dimensions, net height in inches

1 5/16 (33)

External dimensions, net width in inches

23 1/4 (590)

External dimensions, net depth in inches

30 9/16 (775)

External dimensions, gross height in inches

3 3/4 (95)

External dimensions, gross width in inches

24 1/4 (615)

External dimensions, gross depth in inches

30 3/4 (780)

Net weight in lbs

5 (2.13)

Gross weight in lbs

6 (2.90)