


A 605

Module

For optimum reprocessing of vials and test tubes, etc.



- Equipped with 98 injector nozzles
- Max. H 6 5/16" (160 mm), max. Ø 20 mm
- For use in A 503 mobile unit

EAN: 4002516854074 / Material number: 12669870 / Old Material Number: 69660501D

Product category

Large-chamber laboratory glasswashers

•

Machine affiliation

PLW 8615

•

PLW 8616

•

PLW 8617

•

Product category

Insert/module for test tubes

•

Application area

Reprocessing of pipettes/vials

•

Reprocessing of test tubes

•

Reprocessing of laboratory glassware

•

Product characteristics

Material

Stainless steel

Color

Stainless steel

Capacity

Test tubes [number]

98

Phials [number]

98

Dimensions and weight

External dimensions, net height in inches

6 1/4 (158)

External dimensions, net width in inches

10 (254)

External dimensions, net depth in inches

23 13/16 (604)

External dimensions, gross height in inches

7 1/2 (190)

External dimensions, gross width in inches

11 3/8 (288)

External dimensions, gross depth in inches

24 7/16 (620)

Net weight in lbs

6 (2.85)

Gross weight in lbs

7 (3.31)