

A 622

Injector module

For the optimum reprocessing of laboratory glassware, 36 injector positions.



- Can be equipped with up to 36 nozzles based on individual needs
- For use in A 503 and APLW 511 mobile units
- Part of the SimpleLoad system

EAN: 4002516854531 / Material number: 12678490 / Old Material Number: 69662201D

## Product category

Large-chamber laboratory glasswashers

•

## Machine affiliation

PLW 8615

•

PLW 8616

•

PLW 8617

•

## Product category

Insert/module for narrow-neck glassware

•

Insert for various bottles

•

## Application area

Reprocessing of laboratory glassware

•

## Product characteristics

Material

Stainless steel

Color

Stainless steel

Part of the SimpleLoad system

•

## Capacity

Narrow-neck glassware [number]

36

Flasks [number]

36

Round flasks [number]

36

## Dimensions and weight

External dimensions, net height in inches

4 1/2 (113)

External dimensions, net width in inches

8 1/16 (204)

External dimensions, net depth in inches

23 15/16 (608)

External dimensions, gross height in inches

5 3/4 (145)

External dimensions, gross width in inches

8 3/4 (222)

External dimensions, gross depth in inches

24 13/16 (629)

Net weight in lbs

5 (2.23)

Gross weight in lbs

6 (2.61)