

## ProCare Dent 10 MA - 5 l

Liquid detergent, mildly alkaline, 5 l

for the automatic reprocessing of dental and transmission instruments.



- Specifically designed for cleaning dental instruments
- Adapted to the process sequences in Miele washer-disinfectors
- Very gentle on material
- Contains enzymes
- Check for suitability on anodized aluminum prior to use

EAN: 4029577158200 / Material number: 12680720 /  
Old Material Number: 66913204V6

<b>Product category</b>	
Thermal disinfectors	•
<b>Applications</b>	
Removal of blood	•
Removal of organic residues	•
<b>Dispensing</b>	
pH value at 68°F (20°C)	10,1
<b>Dispensing recommendation</b>	
Dispensing recommendation per rinse	3–7 mL/L (0.3–0.7%)
Dispensing recommendation (per rinse)	3–7 mL/L (0.3–0.7%)
Large opening for siphon	•
Intensity of soiling	Light to moderate
Suitable for all levels of water hardness	•
<b>Main applications</b>	
Dental practices	•
<b>Product characteristics</b>	
Medium type	Mildly alkaline Cleaning agent
Aggregate form	Liquid concentrate
Ingredients	Enzymes, Surfactants
Minimum storage temperature in °F/°C	0
Maximum storage temperature in °F/°C	25
<b>Dimensions and weight</b>	
External dimensions, net height in inches	217
External dimensions, net width in inches	142
External dimensions, net depth in inches	248
External dimensions, gross height in inches (mm)	185.00
External dimensions, gross width in inches (mm)	175.00
External dimensions, gross depth in inches (mm)	250.00
Gross weight in lbs (kg)	5.75
Fill volume in ml	5000
<b>Product characteristics</b>	
Medium type	Mildly alkaline Cleaning agent
Aggregate form	Liquid concentrate
Ingredients	Enzymes, Surfactants
Minimum storage temperature in °F/°C	0
Maximum storage temperature in °F/°C	25
Suitable for automatic reprocessing of instruments and utensils	•
Suitable for sensitive materials	•
Suitable for standard instruments and transmission instruments	•
Check for suitability on anodized aluminum prior to use	•
Very gentle on materials	•
<b>Dimensions and weight</b>	
External dimensions, net height in inches	8 9/16 (217)
External dimensions, net width in inches	5 5/8 (142)
External dimensions, net depth in inches	9 13/16 (248)

ProCare Dent 10 MA - 5 l  
Liquid detergent, mildly alkaline, 5 l  
for the automatic reprocessing of dental and transmission instruments.



EAN: 4029577158200 / Material number: 12680720 / Old Material Number: 66913204V6	
External dimensions, gross height in inches	7 5/16 (185)
External dimensions, gross width in inches	6 15/16 (175)
External dimensions, gross depth in inches	9 7/8 (250)
Gross weight in lbs	0.21 (6)
Fill volume in ml	169.07 (5000) fl oz (US) (Milliliter)
Compliance, test certificates, and symbols	
CE	•
GHS/CLP compliant	•